

Calf Implants

by EUROSILICONE®

Ürüne Genel Bakış

- 2 farklı konfigürasyon.
- Simetrik / asimetrik form.
- Yüksek kohesivite.
- Pürüzsüz yüzey.
- İmplantın yerleştirilmesini kolaylaştıran disk.



ES 54


Symmetric Shape



Product Code	Length (cm)	Width (cm)	Projection (cm)	Volume (cc)
54060	16.8	3.30	2.0	60
54090	17.0	3.70	2.8	90
54110	17.7	5.80	2.1	110
54120	18.9	5.20	2.6	120
54160	19.5	6.90	2.3	160
54180	19.0	5.90	3.1	180
54220	18.7	7.10	3.2	220

ES 54

Asymmetric Shape



Product Code	Length (cm)	Width (cm)	Projection (cm)	Volume (cc)
54135	22.0	4.90	2.3	135
54140	18.6	5.30	2.6	140
54150	13.9	6.40	2.9	150
54165	23.9	5.60	2.3	165